



### Chair Backrest Operation

Movement of the Domain chair is accomplished by sliding the black disc/actuator on the foot control. The actuator will slide in four directions; up, down, left and right.

Moving the actuator up or down will raise or lower the chair seat. Moving the actuator left or right will recline or raise the backrest.

The actuator also works in conjunction with the **AP** (Auto Position) button to program the chair. Each of the four directions of actuator movement can represent a preset position for the operator's use. Please refer to the programming instructions.




### Auto Return (Exit Position)

To lower the chair completely and bring the backrest upright, locate the red button labeled **AR**. Press and release this red **AR** button and the chair will move to the exit position.

### Auto Position Selection

To automatically raise move the chair to a pre-set position, press and release the green **AP** button. Immediately slide and then release the black actuator in one of its four directions.



### Programming Auto Position

Move the chair to the desired position using the black actuator. Push and hold the **AP** button. While pressing the **AP** button, slide and hold the black actuator in one of the 4 preset positions (see diagram) until you hear 2 beeps (approx. 5 sec.).

### Programming Last Position

Move the chair to the desired position using the black actuator. Push and hold the **AP** button also push and hold the **LP** button until you hear two beeps.

## Patient Entry or Exit

The armrests may be rotated out of the way for patient entry or exit. Pull up on the release lever under the armrest and rotate the front of the armrest down. To return the armrest swing it back up until you hear it lock in place.

## Chair Rotation

To operate the chair rotation, locate the handle directly behind the seat cushion. To unlock the chair, move the handle clockwise. The chair will swivel 30 degrees left or right for a total of 60 degrees. To lock the chair move the handle counter clockwise.

## Headrest

To adjust the location of the headrest, grasp both sides of the headrest backboard and pull up or push down. The pillow may be pushed down the backboard onto the backrest to accommodate patient requirements.

## Power Light

The green light mounted low on the front of the base shroud indicates that power is reaching the chair controller. To stop the power you may unplug the chair's power cord from the receptacle it is plugged into. If the chair is plugged in and the light does not come on please call your service representative. The chair does have two replaceable fuses rated at 8 amps, 250 volts. These fuses should never fail under normal operation. If a fuse does fail it indicates a problem that should be repaired by a service technician.

## Duty Cycle – Running Time

The chair is intended for the typical intermittent operation of a dental office. Do not run the chair continuously for more than one minute. Also, do not cycle the chair in a manner that would require the chair to move more than 20% of the time.

## Cleaning Instructions

See the included cleaning instructions for the **Boltaflex** or **NaugaSoft** brand upholstery that covers your chair.

## Underwriters Laboratories



DENTAL CHAIR, DOMAIN  
WITH RESPECT TO ELECTRICAL SHOCK, FIRE, MECHANICAL AND  
OTHER SPECIFIED HAZARDS IN ACCORDANCE WITH UL 2601-1 AND  
CAN/CSA C22.2 NO.601.1  
(37RJ)